



GENOSS Spinal Devices

GENOSS 3d Cage™
GENOSS Lumbar Pedicle Screw™






Table of contents

GENOSS 3d Cage™	4
GENOSS Lumbar Pedicle Screw™	8
Order Information	11



GENOSS

3d CageTM


- Excellent bone fusion due to titanium with 3D porous structure
- No subsidence due to similar elastic modulus to cancellous bone
- Less metal artifact by diagnostic imaging

Excellent bone fusion


Osteoconductive porous structure

The unique pore structures and roughened nano-textured surface may help facilitate osseointegration and substantial bone growth throughout the material.


88% porosity



Macro pores



Micro pores




Nano pores

Designed to create better cell attachment and proliferation

GENOSS 3d Cage™ consists of about 88% porosity and average pore size of 800µm (250µm~1,200µm). The porous and micro structure is designed to create a favorable environment for cell attachment and proliferation.

In-vivo test




Beagle dog femur, 12w



3D-Printed porous ti scaffold

Prematured bone



No subsidence

No subsidence due to similar elastic modulus to cancellous bone

GENOSS 3d Cage™ demonstrated better resistance to subsidence than other different materials.

Stress vs. strain curves




Figure 1 Comparison of stress-strain curves of materials common in spinal reconstruction

Elastic modulus under compressive forces






Figure 2 Elastic modulus of all materials


Robert F. Heary, Naresh Parvathreddy, Sujitha Sampath, Nitin Agarwal. Elastic Modulus in the selection of interbody implants[master's thesis].New Jersey(USA): Journal of Spinal Surgery.2017.165p.

Clear visibility

Less metal artifact by diagnostic imaging

High implant porosity provides excellent diagnostic imaging. The fusion area is clearly visible by X-Ray, CT and MRI.






GENOSS

Lumbar Pedicle Screw™

- Enhanced fixation using S.L.A. surface treatment
- Higher bone-to-implant contact
- Faster bone formation on the surface

S.L.A. surface

Enhanced bone fixation performance
using S.L.A. surface treatment



Device description


GENOSS Lumbar Pedicle Screw™ is a device designed for usage in the posterior lumbar inter-body fusion procedures. It consists of pedicle screw(poly, mono type), reduction screw(poly, mono type), set screw, rods, crosslink as well as connecting and bolting elements. Various forms and sizes of these implants are available, so that adaptations can always be made to take into account the pathology and the individual patient.

Indication


GENOSS Lumbar Pedicle Screw™ is intended to provide immobilization and stabilization of spinal segments in the treatment of the following; Degenerative spondylolisthesis objective evidence of neurologic impairment, fracture, dislocation, scoliosis, kyphosis, spinal tumor, and failed previous fusion(pseudarthrosis). In addition, the system is indicated for spinal pedicle screw fixation in patients with severe spondylolisthesis (Grade 3 and 4) at the L5-S1 joint, receiving fusion using autogenous graft, having the device removed after the development of a solid fusion mass.

GENOSS 3d Cage™
Order Information


3d Cage-PLIF Grip sterile packaged cages




Reference No.	Model No.	H(mm)	W(mm)	L(mm)	Θ°
11010001	PL0709220	7	9	22	0
11010002	PL0809220	8	9	22	0
11010003	PL0909220	9	9	22	0
11010004	PL1009220	10	9	22	0
11010005	PL1109220	11	9	22	0
11010006	PL1209220	12	9	22	0
11010007	PL1309220	13	9	22	0
11010008	PL1409220	14	9	22	0
11010009	PL1509220	15	9	22	0
11010010	PL0711220	7	11	22	0
11010011	PL0811220	8	11	22	0
11010012	PL0911220	9	11	22	0
11010013	PL1011220	10	11	22	0
11010014	PL1111220	11	11	22	0
11010015	PL1211220	12	11	22	0
11010016	PL1311220	13	11	22	0
11010017	PL1411220	14	11	22	0
11010018	PL1511220	15	11	22	0
11010019	PL0709224	7	9	22	4
11010020	PL0809224	8	9	22	4
11010021	PL0909224	9	9	22	4
11010022	PL1009224	10	9	22	4
11010023	PL1109224	11	9	22	4
11010024	PL1209224	12	9	22	4
11010025	PL1309224	13	9	22	4
11010026	PL1409224	14	9	22	4
11010027	PL1509224	15	9	22	4
11010028	PL0711224	7	11	22	4
11010029	PL0811224	8	11	22	4
11010030	PL0911224	9	11	22	4
11010031	PL1011224	10	11	22	4
11010032	PL1111224	11	11	22	4
11010033	PL1211224	12	11	22	4
11010034	PL1311224	13	11	22	4
11010035	PL1411224	14	11	22	4
11010036	PL1511224	15	11	22	4
11010037	PL0709228	7	9	22	8
11010038	PL0809228	8	9	22	8
11010039	PL0909228	9	9	22	8
11010040	PL1009228	10	9	22	8
11010041	PL1109228	11	9	22	8
11010042	PL1209228	12	9	22	8
11010043	PL1309228	13	9	22	8
11010044	PL1409228	14	9	22	8
11010045	PL1509228	15	9	22	8
11010046	PL0711228	7	11	22	8
11010047	PL0811228	8	11	22	8
11010048	PL0911228	9	11	22	8
11010049	PL1011228	10	11	22	8
11010050	PL1111228	11	11	22	8
11010051	PL1211228	12	11	22	8
11010052	PL1311228	13	11	22	8
11010053	PL1411228	14	11	22	8
11010054	PL1511228	15	11	22	8




Reference No.	Model No.	H(mm)	W(mm)	L(mm)	Θ°
11010055	PL0709240	7	9	24	0
11010056	PL0809240	8	9	24	0
11010057	PL0909240	9	9	24	0
11010058	PL1009240	10	9	24	0
11010059	PL1109240	11	9	24	0
11010060	PL1209240	12	9	24	0
11010061	PL1309240	13	9	24	0
11010062	PL1409240	14	9	24	0
11010063	PL1509240	15	9	24	0
11010064	PL0711240	7	11	24	0
11010065	PL0811240	8	11	24	0
11010066	PL0911240	9	11	24	0
11010067	PL1011240	10	11	24	0
11010068	PL1111240	11	11	24	0
11010069	PL1211240	12	11	24	0
11010070	PL1311240	13	11	24	0
11010071	PL1411240	14	11	24	0
11010072	PL1511240	15	11	24	0
11010073	PL0709244	7	9	24	4
11010074	PL0809244	8	9	24	4
11010075	PL0909244	9	9	24	4
11010076	PL1009244	10	9	24	4
11010077	PL1109244	11	9	24	4
11010078	PL1209244	12	9	24	4
11010079	PL1309244	13	9	24	4
11010080	PL1409244	14	9	24	4
11010081	PL1509244	15	9	24	4
11010082	PL0711244	7	11	24	4
11010083	PL0811244	8	11	24	4
11010084	PL0911244	9	11	24	4
11010085	PL1011244	10	11	24	4
11010086	PL1111244	11	11	24	4
11010087	PL1211244	12	11	24	4
11010088	PL1311244	13	11	24	4
11010089	PL1411244	14	11	24	4
11010090	PL1511244	15	11	24	4
11010091	PL0709248	7	9	24	8
11010092	PL0809248	8	9	24	8
11010093	PL0909248	9	9	24	8
11010094	PL1009248	10	9	24	8
11010095	PL1109248	11	9	24	8
11010096	PL1209248	12	9	24	8
11010097	PL1309248	13	9	24	8
11010098	PL1409248	14	9	24	8
11010099	PL1509248	15	9	24	8
11010100	PL0711248	7	11	24	8
11010101	PL0811248	8	11	24	8
11010102	PL0911248	9	11	24	8
11010103	PL1011248	10	11	24	8
11010104	PL1111248	11	11	24	8
11010105	PL1211248	12	11	24	8
11010106	PL1311248	13	11	24	8
11010107	PL1411248	14	11	24	8
11010108	PL1511248	15	11	24	8




Reference No.	Model No.	H(mm)	W(mm)	L(mm)	Θ°
11010109	PL0709260	7	9	26	0
11010110	PL0809260	8	9	26	0
11010111	PL0909260	9	9	26	0
11010112	PL1009260	10	9	26	0
11010113	PL1109260	11	9	26	0
11010114	PL1209260	12	9	26	0
11010115	PL1309260	13	9	26	0
11010116	PL1409260	14	9	26	0
11010117	PL1509260	15	9	26	0
11010118	PL0711260	7	11	26	0
11010119	PL0811260	8	11	26	0
11010120	PL0911260	9	11	26	0
11010121	PL1011260	10	11	26	0
11010122	PL1111260	11	11	26	0
11010123	PL1211260	12	11	26	0
11010124	PL1311260	13	11	26	0
11010125	PL1411260	14	11	26	0
11010126	PL1511260	15	11	26	0
11010127	PL0709264	7	9	26	4
11010128	PL0809264	8	9	26	4
11010129	PL0909264	9	9	26	4
11010130	PL1009264	10	9	26	4
11010131	PL1109264	11	9	26	4
11010132	PL1209264	12	9	26	4
11010133	PL1309264	13	9	26	4
11010134	PL1409264	14	9	26	4
11010135	PL1509264	15	9	26	4
11010136	PL0711264	7	11	26	4
11010137	PL0811264	8	11	26	4
11010138	PL0911264	9	11	26	4
11010139	PL1011264	10	11	26	4
11010140	PL1111264	11	11	26	4
11010141	PL1211264	12	11	26	4
11010142	PL1311264	13	11	26	4
11010143	PL1411264	14	11	26	4
11010144	PL1511264	15	11	26	4
11010145	PL0709268	7	9	26	8
11010146	PL0809268	8	9	26	8
11010147	PL0909268	9	9	26	8
11010148	PL1009268	10	9	26	8
11010149	PL1109268	11	9	26	8
11010150	PL1209268	12	9	26	8
11010151	PL1309268	13	9	26	8
11010152	PL1409268	14	9	26	8
11010153	PL1509268	15	9	26	8
11010154	PL0711268	7	11	26	8
11010155	PL0811268	8	11	26	8
11010156	PL0911268	9	11	26	8
11010157	PL1011268	10	11	26	8
11010158	PL1111268	11	11	26	8
11010159	PL1211268	12	11	26	8
11010160	PL1311268	13	11	26	8
11010161	PL1411268	14	11	26	8
11010162	PL1511268	15	11	26	8



Reference No.	Model No.	H(mm)	W(mm)	L(mm)	Θ°
11010163	PL0709280	7	9	28	0
11010164	PL0809280	8	9	28	0
11010165	PL0909280	9	9	28	0
11010166	PL1009280	10	9	28	0
11010167	PL1109280	11	9	28	0
11010168	PL1209280	12	9	28	0
11010169	PL1309280	13	9	28	0
11010170	PL1409280	14	9	28	0
11010171	PL1509280	15	9	28	0
11010172	PL0711280	7	11	28	0
11010173	PL0811280	8	11	28	0
11010174	PL0911280	9	11	28	0
11010175	PL1011280	10	11	28	0
11010176	PL1111280	11	11	28	0
11010177	PL1211280	12	11	28	0
11010178	PL1311280	13	11	28	0
11010179	PL1411280	14	11	28	0
11010180	PL1511280	15	11	28	0
11010181	PL0709284	7	9	28	4
11010182	PL0809284	8	9	28	4
11010183	PL0909284	9	9	28	4
11010184	PL1009284	10	9	28	4
11010185	PL1109284	11	9	28	4
11010186	PL1209284	12	9	28	4
11010187	PL1309284	13	9	28	4
11010188	PL1409284	14	9	28	4
11010189	PL1509284	15	9	28	4
11010190	PL0711284	7	11	28	4
11010191	PL0811284	8	11	28	4
11010192	PL0911284	9	11	28	4
11010193	PL1011284	10	11	28	4
11010194	PL1111284	11	11	28	4
11010195	PL1211284	12	11	28	4
11010196	PL1311284	13	11	28	4
11010197	PL1411284	14	11	28	4
11010198	PL1511284	15	11	28	4
11010199	PL0709288	7	9	28	8
11010200	PL0809288	8	9	28	8
11010201	PL0909288	9	9	28	8
11010202	PL1009288	10	9	28	8
11010203	PL1109288	11	9	28	8
11010204	PL1209288	12	9	28	8
11010205	PL1309288	13	9	28	8
11010206	PL1409288	14	9	28	8
11010207	PL1509288	15	9	28	8
11010208	PL0711288	7	11	28	8
11010209	PL0811288	8	11	28	8
11010210	PL0911288	9	11	28	8
11010211	PL1011288	10	11	28	8
11010212	PL1111288	11	11	28	8
11010213	PL1211288	12	11	28	8
11010214	PL1311288	13	11	28	8
11010215	PL1411288	14	11	28	8
11010216	PL1511288	15	11	28	8




Reference No.	Model No.	H(mm)	W(mm)	L(mm)	Θ°
11010217	PL0709300	7	9	30	0
11010218	PL0809300	8	9	30	0
11010219	PL0909300	9	9	30	0
11010220	PL1009300	10	9	30	0
11010221	PL1109300	11	9	30	0
11010222	PL1209300	12	9	30	0
11010223	PL1309300	13	9	30	0
11010224	PL1409300	14	9	30	0
11010225	PL1509300	15	9	30	0
11010226	PL0711300	7	11	30	0
11010227	PL0811300	8	11	30	0
11010228	PL0911300	9	11	30	0
11010229	PL1011300	10	11	30	0
11010230	PL1111300	11	11	30	0
11010231	PL1211300	12	11	30	0
11010232	PL1311300	13	11	30	0
11010233	PL1411300	14	11	30	0
11010234	PL1511300	15	11	30	0
11010235	PL0709304	7	9	30	4
11010236	PL0809304	8	9	30	4
11010237	PL0909304	9	9	30	4
11010238	PL1009304	10	9	30	4
11010239	PL1109304	11	9	30	4
11010240	PL1209304	12	9	30	4
11010241	PL1309304	13	9	30	4
11010242	PL1409304	14	9	30	4
11010243	PL1509304	15	9	30	4
11010244	PL0711304	7	11	30	4
11010245	PL0811304	8	11	30	4
11010246	PL0911304	9	11	30	4
11010247	PL1011304	10	11	30	4
11010248	PL1111304	11	11	30	4
11010249	PL1211304	12	11	30	4
11010250	PL1311304	13	11	30	4
11010251	PL1411304	14	11	30	4
11010252	PL1511304	15	11	30	4
11010253	PL0709308	7	9	30	8
11010254	PL0809308	8	9	30	8
11010255	PL0909308	9	9	30	8
11010256	PL1009308	10	9	30	8
11010257	PL1109308	11	9	30	8
11010258	PL1209308	12	9	30	8
11010259	PL1309308	13	9	30	8
11010260	PL1409308	14	9	30	8
11010261	PL1509308	15	9	30	8
11010262	PL0711308	7	11	30	8
11010263	PL0811308	8	11	30	8
11010264	PL0911308	9	11	30	8
11010265	PL1011308	10	11	30	8
11010266	PL1111308	11	11	30	8
11010267	PL1211308	12	11	30	8
11010268	PL1311308	13	11	30	8
11010269	PL1411308	14	11	30	8
11010270	PL1511308	15	11	30	8




Reference No.	Model No.	H(mm)	W(mm)	L(mm)	Θ°
11010271	PL0709320	7	9	32	0
11010272	PL0809320	8	9	32	0
11010273	PL0909320	9	9	32	0
11010274	PL1009320	10	9	32	0
11010275	PL1109320	11	9	32	0
11010276	PL1209320	12	9	32	0
11010277	PL1309320	13	9	32	0
11010278	PL1409320	14	9	32	0
11010279	PL1509320	15	9	32	0
11010280	PL0711320	7	11	32	0
11010281	PL0811320	8	11	32	0
11010282	PL0911320	9	11	32	0
11010283	PL1011320	10	11	32	0
11010284	PL1111320	11	11	32	0
11010285	PL1211320	12	11	32	0
11010286	PL1311320	13	11	32	0
11010287	PL1411320	14	11	32	0
11010288	PL1511320	15	11	32	0
11010289	PL0709324	7	9	32	4
11010290	PL0809324	8	9	32	4
11010291	PL0909324	9	9	32	4
11010292	PL1009324	10	9	32	4
11010293	PL1109324	11	9	32	4
11010294	PL1209324	12	9	32	4
11010295	PL1309324	13	9	32	4
11010296	PL1409324	14	9	32	4
11010297	PL1509324	15	9	32	4
11010298	PL0711324	7	11	32	4
11010299	PL0811324	8	11	32	4
11010300	PL0911324	9	11	32	4
11010301	PL1011324	10	11	32	4
11010302	PL1111324	11	11	32	4
11010303	PL1211324	12	11	32	4
11010304	PL1311324	13	11	32	4
11010305	PL1411324	14	11	32	4
11010306	PL1511324	15	11	32	4
11010307	PL0709328	7	9	32	8
11010308	PL0809328	8	9	32	8
11010309	PL0909328	9	9	32	8
11010310	PL1009328	10	9	32	8
11010311	PL1109328	11	9	32	8
11010312	PL1209328	12	9	32	8
11010313	PL1309328	13	9	32	8
11010314	PL1409328	14	9	32	8
11010315	PL1509328	15	9	32	8
11010316	PL0711328	7	11	32	8
11010317	PL0811328	8	11	32	8
11010318	PL0911328	9	11	32	8
11010319	PL1011328	10	11	32	8
11010320	PL1111328	11	11	32	8
11010321	PL1211328	12	11	32	8
11010322	PL1311328	13	11	32	8
11010323	PL1411328	14	11	32	8
11010324	PL1511328	15	11	32	8


3d Cage-TLIF Grip sterile packaged cages




Reference No.	Model No.	H(mm)	W(mm)	L(mm)	θ°
11020001	TL0710280	7	10	28	0
11020002	TL0810280	8	10	28	0
11020003	TL0910280	9	10	28	0
11020004	TL1010280	10	10	28	0
11020005	TL1110280	11	10	28	0
11020006	TL1210280	12	10	28	0
11020007	TL1310280	13	10	28	0
11020008	TL1410280	14	10	28	0
11020009	TL1510280	15	10	28	0
11020010	TL0712280	7	12	28	0
11020011	TL0812280	8	12	28	0
11020012	TL0912280	9	12	28	0
11020013	TL1012280	10	12	28	0
11020014	TL1112280	11	12	28	0
11020015	TL1212280	12	12	28	0
11020016	TL1312280	13	12	28	0
11020017	TL1412280	14	12	28	0
11020018	TL1512280	15	12	28	0
11020019	TL0710284	7	10	28	4
11020020	TL0810284	8	10	28	4
11020021	TL0910284	9	10	28	4
11020022	TL1010284	10	10	28	4
11020023	TL1110284	11	10	28	4
11020024	TL1210284	12	10	28	4
11020025	TL1310284	13	10	28	4
11020026	TL1410284	14	10	28	4
11020027	TL1510284	15	10	28	4
11020028	TL0712284	7	12	28	4
11020029	TL0812284	8	12	28	4
11020030	TL0912284	9	12	28	4
11020031	TL1012284	10	12	28	4
11020032	TL1112284	11	12	28	4
11020033	TL1212284	12	12	28	4
11020034	TL1312284	13	12	28	4
11020035	TL1412284	14	12	28	4
11020036	TL1512284	15	12	28	4
11020037	TL0710288	7	10	28	8
11020038	TL0810288	8	10	28	8
11020039	TL0910288	9	10	28	8
11020040	TL1010288	10	10	28	8
11020041	TL1110288	11	10	28	8
11020042	TL1210288	12	10	28	8
11020043	TL1310288	13	10	28	8
11020044	TL1410288	14	10	28	8
11020045	TL1510288	15	10	28	8
11020046	TL0712288	7	12	28	8
11020047	TL0812288	8	12	28	8
11020048	TL0912288	9	12	28	8
11020049	TL1012288	10	12	28	8
11020050	TL1112288	11	12	28	8
11020051	TL1212288	12	12	28	8
11020052	TL1312288	13	12	28	8
11020053	TL1412288	14	12	28	8
11020054	TL1512288	15	12	28	8
11020055	TL0710282	7	10	28	12
11020056	TL0810282	8	10	28	12
11020057	TL0910282	9	10	28	12
11020058	TL1010282	10	10	28	12
11020059	TL1110282	11	10	28	12
11020060	TL1210282	12	10	28	12
11020061	TL1310282	13	10	28	12
11020062	TL1410282	14	10	28	12
11020063	TL1510282	15	10	28	12
11020064	TL0712282	7	12	28	12
11020065	TL0812282	8	12	28	12
11020066	TL0912282	9	12	28	12
11020067	TL1012282	10	12	28	12
11020068	TL1112282	11	12	28	12
11020069	TL1212282	12	12	28	12
11020070	TL1312282	13	12	28	12
11020071	TL1412282	14	12	28	12
11020072	TL1512282	15	12	28	12



Reference No.	Model No.	H(mm)	W(mm)	L(mm)	θ°
11020073	TL0710300	7	10	30	0
11020074	TL0810300	8	10	30	0
11020075	TL0910300	9	10	30	0
11020076	TL1010300	10	10	30	0
11020077	TL1110300	11	10	30	0
11020078	TL1210300	12	10	30	0
11020079	TL1310300	13	10	30	0
11020080	TL1410300	14	10	30	0
11020081	TL1510300	15	10	30	0
11020082	TL0712300	7	12	30	0
11020083	TL0812300	8	12	30	0
11020084	TL0912300	9	12	30	0
11020085	TL1012300	10	12	30	0
11020086	TL1112300	11	12	30	0
11020087	TL1212300	12	12	30	0
11020088	TL1312300	13	12	30	0
11020089	TL1412300	14	12	30	0
11020090	TL1512300	15	12	30	0
11020091	TL0710304	7	10	30	4
11020092	TL0810304	8	10	30	4
11020093	TL0910304	9	10	30	4
11020094	TL1010304	10	10	30	4
11020095	TL1110304	11	10	30	4
11020096	TL1210304	12	10	30	4
11020097	TL1310304	13	10	30	4
11020098	TL1410304	14	10	30	4
11020099	TL1510304	15	10	30	4
11020100	TL0712304	7	12	30	4
11020101	TL0812304	8	12	30	4
11020102	TL0912304	9	12	30	4
11020103	TL1012304	10	12	30	4
11020104	TL1112304	11	12	30	4
11020105	TL1212304	12	12	30	4
11020106	TL1312304	13	12	30	4
11020107	TL1412304	14	12	30	4
11020108	TL1512304	15	12	30	4
11020109	TL0710308	7	10	30	8
11020110	TL0810308	8	10	30	8
11020111	TL0910308	9	10	30	8
11020112	TL1010308	10	10	30	8
11020113	TL1110308	11	10	30	8
11020114	TL1210308	12	10	30	8
11020115	TL1310308	13	10	30	8
11020116	TL1410308	14	10	30	8
11020117	TL1510308	15	10	30	8
11020118	TL0712308	7	12	30	8
11020119	TL0812308	8	12	30	8
11020120	TL0912308	9	12	30	8
11020121	TL1012308	10	12	30	8
11020122	TL1112308	11	12	30	8
11020123	TL1212308	12	12	30	8
11020124	TL1312308	13	12	30	8
11020125	TL1412308	14	12	30	8
11020126	TL1512308	15	12	30	8
11020127	TL0710302	7	10	30	12
11020128	TL0810302	8	10	30	12
11020129	TL0910302	9	10	30	12
11020130	TL1010302	10	10	30	12
11020131	TL1110302	11	10	30	12
11020132	TL1210302	12	10	30	12
11020133	TL1310302	13	10	30	12
11020134	TL1410302	14	10	30	12
11020135	TL1510302	15	10	30	12
11020136	TL0712302	7	12	30	12
11020137	TL0812302	8	12	30	12
11020138	TL0912302	9	12	30	12
11020139	TL1012302	10	12	30	12
11020140	TL1112302	11	12	30	12
11020141	TL1212302	12	12	30	12
11020142	TL1312302	13	12	30	12
11020143	TL1412302	14	12	30	12
11020144	TL1512302	15	12	30	12




Reference No.	Model No.	H(mm)	W(mm)	L(mm)	θ°
11020145	TL0710320	7	10	32	0
11020146	TL0810320	8	10	32	0
11020147	TL0910320	9	10	32	0
11020148	TL1010320	10	10	32	0
11020149	TL1110320	11	10	32	0
11020150	TL1210320	12	10	32	0
11020151	TL1310320	13	10	32	0
11020152	TL1410320	14	10	32	0
11020153	TL1510320	15	10	32	0
11020154	TL0712320	7	12	32	0
11020155	TL0812320	8	12	32	0
11020156	TL0912320	9	12	32	0
11020157	TL1012320	10	12	32	0
11020158	TL1112320	11	12	32	0
11020159	TL1212320	12	12	32	0
11020160	TL1312320	13	12	32	0
11020161	TL1412320	14	12	32	0
11020162	TL1512320	15	12	32	0
11020163	TL0710324	7	10	32	4
11020164	TL0810324	8	10	32	4
11020165	TL0910324	9	10	32	4
11020166	TL1010324	10	10	32	4
11020167	TL1110324	11	10	32	4
11020168	TL1210324	12	10	32	4
11020169	TL1310324	13	10	32	4
11020170	TL1410324	14	10	32	4
11020171	TL1510324	15	10	32	4
11020172	TL0712324	7	12	32	4
11020173	TL0812324	8	12	32	4
11020174	TL0912324	9	12	32	4
11020175	TL1012324	10	12	32	4
11020176	TL1112324	11	12	32	4
11020177	TL1212324	12	12	32	4
11020178	TL1312324	13	12	32	4
11020179	TL1412324	14	12	32	4
11020180	TL1512324	15	12	32	4
11020181	TL0710328	7	10	32	8
11020182	TL0810328	8	10	32	8
11020183	TL0910328	9	10	32	8
11020184	TL1010328	10	10	32	8
11020185	TL1110328	11	10	32	8
11020186	TL1210328	12	10	32	8
11020187	TL1310328	13	10	32	8
11020188	TL1410328	14	10	32	8
11020189	TL1510328	15	10	32	8
11020190	TL0712328	7	12	32	8
11020191	TL0812328	8	12	32	8
11020192	TL0912328	9	12	32	8
11020193	TL1012328	10	12	32	8
11020194	TL1112328	11	12	32	8
11020195	TL1212328	12	12	32	8
11020196	TL1312328	13	12	32	8
11020197	TL1412328	14	12	32	8
11020198	TL1512328	15	12	32	8
11020199	TL0710322	7	10	32	12
11020200	TL0810322	8	10	32	12
11020201	TL0910322	9	10	32	12
11020202	TL1010322	10	10	32	12
11020203	TL1110322	11	10	32	12
11020204	TL1210322	12	10	32	12
11020205	TL1310322	13	10	32	12
11020206	TL1410322	14	10	32	12
11020207	TL1510322	15	10	32	12
11020208	TL0712322	7	12	32	12
11020209	TL0812322	8	12	32	12
11020210	TL0912322	9	12	32	12
11020211	TL1012322	10	12	32	12
11020212	TL1112322	11	12	32	12
11020213	TL1212322	12	12	32	12
11020214	TL1312322	13	12	32	12
11020215	TL1412322	14	12	32	12
11020216	TL1512322	15	12	32	12




Reference No.


Model No.	H(mm)	W(mm)	L(mm)	θ°	
11020217	TL0710340	7	10	34	0
11020218	TL0810340	8	10	34	0
11020219	TL0910340	9	10	34	0
11020220	TL1010340	10	10	34	0
11020221	TL1110340	11	10	34	0
11020222	TL1210340	12	10	34	0
11020223	TL1310340	13	10	34	0
11020224	TL1410340	14	10	34	0
11020225	TL1510340	15	10	34	0
11020226	TL0712340	7	12	34	0
11020227	TL0812340	8	12	34	0
11020228	TL0912340	9	12	34	0
11020229	TL1012340	10	12	34	0
11020230	TL1112340	11	12	34	0
11020231	TL1212340	12	12	34	0
11020232	TL1312340	13	12	34	0
11020233	TL1412340	14	12	34	0
11020234	TL1512340	15	12	34	0
11020235	TL0710344	7	10	34	4
11020236	TL0810344	8	10	34	4
11020237	TL0910344	9	10	34	4
11020238	TL1010344	10	10	34	4
11020239	TL1110344	11	10	34	4
11020240	TL1210344	12	10	34	4
11020241	TL1310344	13	10	34	4
11020242	TL1410344	14	10	34	4
11020243	TL1510344	15	10	34	4
11020244	TL0712344	7	12	34	4
11020245	TL0812344	8	12	34	4
11020246	TL0912344	9	12	34	4
11020247	TL1012344	10	12	34	4
11020248	TL1112344	11	12	34	4
11020249	TL1212344	12	12	34	4
11020250	TL1312344	13	12	34	4
11020251	TL1412344	14	12	34	4
11020252	TL1512344	15	12	34	4
11020253	TL0710348	7	10	34	8
11020254	TL0810348	8	10	34	8
11020255	TL0910348	9	10	34	8
11020256	TL1010348	10	10	34	8
11020257	TL1110348	11	10	34	8
11020258	TL1210348	12	10	34	8
11020259	TL1310348	13	10	34	8
11020260	TL1410348	14	10	34	8
11020261	TL1510348	15	10	34	8
11020262	TL0712348	7	12	34	8
11020263	TL0812348	8	12	34	8
11020264	TL0912348	9	12	34	8
11020265	TL1012348	10	12	34	8
11020266	TL1112348	11	12	34	8
11020267	TL1212348	12	12	34	8
11020268	TL1312348	13	12	34	8
11020269	TL1412348	14	12	34	8
11020270	TL1512348	15	12	34	8
11020271	TL0710342	7	10	34	12
11020272	TL0810342	8	10	34	12
11020273	TL0910342	9	10	34	12
11020274	TL1010342	10	10	34	12
11020275	TL1110342	11	10	34	12
11020276	TL1210342	12	10	34	12
11020277	TL1310342	13	10	34	12
11020278	TL1410342	14	10	34	12
11020279	TL1510342	15	10	34	12
11020280	TL0712342	7	12	34	12
11020281	TL0812342	8	12	34	12
11020282	TL0912342	9	12	34	12
11020283	TL1012342	10	12	34	12
11020284	TL1112342	11	12	34	12
11020285	TL1212342	12	12	34	12
11020286	TL1312342	13	12	34	12
11020287	TL1412342	14	12	34	12
11020288	TL1512342	15	12	34	12




Reference No.	Model No.	H(mm)	W(mm)	L(mm)	θ°
11020289	TL0710360	7	10	36	0
11020290	TL0810360	8	10	36	0
11020291	TL0910360	9	10	36	0
11020292	TL1010360	10	10	36	0
11020293	TL1110360	11	10	36	0
11020294	TL1210360	12	10	36	0
11020295	TL1310360	13	10	36	0
11020296	TL1410360	14	10	36	0
11020297	TL1510360	15	10	36	0
11020298	TL0712360	7	12	36	0
11020299	TL0812360	8	12	36	0
11020300	TL0912360	9	12	36	0
11020301	TL1012360	10	12	36	0
11020302	TL1112360	11	12	36	0
11020303	TL1212360	12	12	36	0
11020304	TL1312360	13	12	36	0
11020305	TL1412360	14	12	36	0
11020306	TL1512360	15	12	36	0
11020307	TL0710364	7	10	36	4
11020308	TL0810364	8	10	36	4
11020309	TL0910364	9	10	36	4
11020310	TL1010364	10	10	36	4
11020311	TL1110364	11	10	36	4
11020312	TL1210364	12	10	36	4
11020313	TL1310364	13	10	36	4
11020314	TL1410364	14	10	36	4
11020315	TL1510364	15	10	36	4
11020316	TL0712364	7	12	36	4
11020317	TL0812364	8	12	36	4
11020318	TL0912364	9	12	36	4
11020319	TL1012364	10	12	36	4
11020320	TL1112364	11	12	36	4
11020321	TL1212364	12	12	36	4
11020322	TL1312364	13	12	36	4
11020323	TL1412364	14	12	36	4
11020324	TL1512364	15	12	36	4
11020325	TL0710368	7	10	36	8
11020326	TL0810368	8	10	36	8
11020327	TL0910368	9	10	36	8
11020328	TL1010368	10	10	36	8
11020329	TL1110368	11	10	36	8
11020330	TL1210368	12	10	36	8
11020331	TL1310368	13	10	36	8
11020332	TL1410368	14	10	36	8
11020333	TL1510368	15	10	36	8
11020334	TL0712368	7	12	36	8
11020335	TL0812368	8	12	36	8
11020336	TL0912368	9	12	36	8
11020337	TL1012368	10	12	36	8
11020338	TL1112368	11	12	36	8
11020339	TL1212368	12	12	36	8
11020340	TL1312368	13	12	36	8
11020341	TL1412368	14	12	36	8
11020342	TL1512368	15	12	36	8
11020343	TL0710362	7	10	36	12
11020344	TL0810362	8	10	36	12
11020345	TL0910362	9	10	36	12
11020346	TL1010362	10	10	36	12
11020347	TL1110362	11	10	36	12
11020348	TL1210362	12	10	36	12
11020349	TL1310362	13	10	36	12
11020350	TL1410362	14	10	36	12
11020351	TL1510362	15	10	36	12
11020352	TL0712362	7	12	36	12
11020353	TL0812362	8	12	36	12
11020354	TL0912362	9	12	36	12
11020355	TL1012362	10	12	36	12
11020356	TL1112362	11	12	36	12
11020357	TL1212362	12	12	36	12
11020358	TL1312362	13	12	36	12
11020359	TL1412362	14	12	36	12
11020360	TL1512362	15	12	36	12




Reference No.	Model No.	H(mm)	W(mm)	L(mm)	θ°
11020361	TL0710380	7	10	38	0
11020362	TL0810380	8	10	38	0
11020363	TL0910380	9	10	38	0
11020364	TL1010380	10	10	38	0
11020365	TL1110380	11	10	38	0
11020366	TL1210380	12	10	38	0
11020367	TL1310380	13	10	38	0
11020368	TL1410380	14	10	38	0
11020369	TL1510380	15	10	38	0
11020370	TL0712380	7	12	38	0
11020371	TL0812380	8	12	38	0
11020372	TL0912380	9	12	38	0
11020373	TL1012380	10	12	38	0
11020374	TL1112380	11	12	38	0
11020375	TL1212380	12	12	38	0
11020376	TL1312380	13	12	38	0
11020377	TL1412380	14	12	38	0
11020378	TL1512380	15	12	38	0
11020379	TL0710384	7	10	38	4
11020380	TL0810384	8	10	38	4
11020381	TL0910384	9	10	38	4
11020382	TL1010384	10	10	38	4
11020383	TL1110384	11	10	38	4
11020384	TL1210384	12	10	38	4
11020385	TL1310384	13	10	38	4
11020386	TL1410384	14	10	38	4
11020387	TL1510384	15	10	38	4
11020388	TL0712384	7	12	38	4
11020389	TL0812384	8	12	38	4
11020390	TL0912384	9	12	38	4
11020391	TL1012384	10	12	38	4
11020392	TL1112384	11	12	38	4
11020393	TL1212384	12	12	38	4
11020394	TL1312384	13	12	38	4
11020395	TL1412384	14	12	38	4
11020396	TL1512384	15	12	38	4
11020397	TL0710388	7	10	38	8
11020398	TL0810388	8	10	38	8
11020399	TL0910388	9	10	38	8
11020400	TL1010388	10	10	38	8
11020401	TL1110388	11	10	38	8
11020402	TL1210388	12	10	38	8
11020403	TL1310388	13	10	38	8
11020404	TL1410388	14	10	38	8
11020405	TL1510388	15	10	38	8
11020406	TL0712388	7	12	38	8
11020407	TL0812388	8	12	38	8
11020408	TL0912388	9	12	38	8
11020409	TL1012388	10	12	38	8
11020410	TL1112388	11	12	38	8
11020411	TL1212388	12	12	38	8
11020412	TL1312388	13	12	38	8
11020413	TL1412388	14	12	38	8
11020414	TL1512388	15	12	38	8
11020415	TL0710382	7	10	38	12
11020416	TL0810382	8	10	38	12
11020417	TL0910382	9	10	38	12
11020418	TL1010382	10	10	38	12
11020419	TL1110382	11	10	38	12
11020420	TL1210382	12	10	38	12
11020421	TL1310382	13	10	38	12
11020422	TL1410382	14	10	38	12
11020423	TL1510382	15	10	38	12
11020424	TL0712382	7	12	38	12
11020425	TL0812382	8	12	38	12
11020426	TL0912382	9	12	38	12
11020427	TL1012382	10	12	38	12
11020428	TL1112382	11	12	38	12
11020429	TL1212382	12	12	38	12
11020430	TL1312382	13	12	38	12
11020431	TL1412382	14	12	38	12
11020432	TL1512382	15	12	38	12




Reference No.	Model No.	H(mm)	W(mm)	L(mm)	θ°
11020433	TL0710400	7	10	40	0
11020434	TL0810400	8	10	40	0
11020435	TL0910400	9	10	40	0
11020436	TL1010400	10	10	40	0
11020437	TL1110400	11	10	40	0
11020438	TL1210400	12	10	40	0
11020439	TL1310400	13	10	40	0
11020440	TL1410400	14	10	40	0
11020441	TL1510400	15	10	40	0
11020442	TL0712400	7	12	40	0
11020443	TL0812400	8	12	40	0
11020444	TL0912400	9	12	40	0
11020445	TL1012400	10	12	40	0
11020446	TL1112400	11	12	40	0
11020447	TL1212400	12	12	40	0
11020448	TL1312400	13	12	40	0
11020449	TL1412400	14	12	40	0
11020450	TL1512400	15	12	40	0
11020451	TL0710404	7	10	40	4
11020452	TL0810404	8	10	40	4
11020453	TL0910404	9	10	40	4
11020454	TL1010404	10	10	40	4
11020455	TL1110404	11	10	40	4
11020456	TL1210404	12	10	40	4
11020457	TL1310404	13	10	40	4
11020458	TL1410404	14	10	40	4
11020459	TL1510404	15	10	40	4
11020460	TL0712404	7	12	40	4
11020461	TL0812404	8	12	40	4
11020462	TL0912404	9	12	40	4
11020463	TL1012404	10	12	40	4
11020464	TL1112404	11	12	40	4
11020465	TL1212404	12	12	40	4
11020466	TL1312404	13	12	40	4
11020467	TL1412404	14	12	40	4
11020468	TL1512404	15	12	40	4
11020469	TL0710408	7	10	40	8
11020470	TL0810408	8	10	40	8
11020471	TL0910408	9	10	40	8
11020472	TL1010408	10	10	40	8
11020473	TL1110408	11	10	40	8
11020474	TL1210408	12	10	40	8
11020475	TL1310408	13	10	40	8
11020476	TL1410408	14	10	40	8
11020477	TL1510408	15	10	40	8
11020478	TL0712408	7	12	40	8
11020479	TL0812408	8	12	40	8
11020480	TL0912408	9	12	40	8
11020481	TL1012408	10	12	40	8
11020482	TL1112408	11	12	40	8
11020483	TL1212408	12	12	40	8
11020484	TL1312408	13	12	40	8
11020485	TL1412408	14	12	40	8
11020486	TL1512408	15	12	40	8
11020487	TL0710402	7	10	40	12
11020488	TL0810402	8	10	40	12
11020489	TL0910402	9	10	40	12
11020490	TL1010402	10	10	40	12
11020491	TL1110402	11	10	40	12
11020492	TL1210402	12	10	40	12
11020493	TL1310402	13	10	40	12
11020494	TL1410402	14	10	40	12
11020495	TL1510402	15	10	40	12
11020496	TL0712402	7	12	40	12
11020497	TL0812402	8	12	40	12
11020498	TL0912402	9	12	40	12
11020499	TL1012402	10	12	40	12
11020500	TL1112402	11	12	40	12
11020501	TL1212402	12	12	40	12
11020502	TL1312402	13	12	40	12
11020503	TL1412402	14	12	40	12
11020504	TL1512402	15	12	40	12




3d Cage-PLIF sterile packaged cages




Reference No.	Model No.	H(mm)	W(mm)	L(mm)	Θ°
11020505	PT0709220	7	9	22	0
11020506	PT0809220	8	9	22	0
11020507	PT0909220	9	9	22	0
11020508	PT1009220	10	9	22	0
11020509	PT1109220	11	9	22	0
11020510	PT1209220	12	9	22	0
11020511	PT1309220	13	9	22	0
11020512	PT1409220	14	9	22	0
11020513	PT1509220	15	9	22	0
11020514	PT0711220	7	11	22	0
11020515	PT0811220	8	11	22	0
11020516	PT0911220	9	11	22	0
11020517	PT1011220	10	11	22	0
11020518	PT1111220	11	11	22	0
11020519	PT1211220	12	11	22	0
11020520	PT1311220	13	11	22	0
11020521	PT1411220	14	11	22	0
11020522	PT1511220	15	11	22	0
11020523	PT0709224	7	9	22	4
11020524	PT0809224	8	9	22	4
11020525	PT0909224	9	9	22	4
11020526	PT1009224	10	9	22	4
11020527	PT1109224	11	9	22	4
11020528	PT1209224	12	9	22	4
11020529	PT1309224	13	9	22	4
11020530	PT1409224	14	9	22	4
11020531	PT1509224	15	9	22	4
11020532	PT0711224	7	11	22	4
11020533	PT0811224	8	11	22	4
11020534	PT0911224	9	11	22	4
11020535	PT1011224	10	11	22	4
11020536	PT1111224	11	11	22	4
11020537	PT1211224	12	11	22	4
11020538	PT1311224	13	11	22	4
11020539	PT1411224	14	11	22	4
11020540	PT1511224	15	11	22	4
11020541	PT0709228	7	9	22	8
11020542	PT0809228	8	9	22	8
11020543	PT0909228	9	9	22	8
11020544	PT1009228	10	9	22	8
11020545	PT1109228	11	9	22	8
11020546	PT1209228	12	9	22	8
11020547	PT1309228	13	9	22	8
11020548	PT1409228	14	9	22	8
11020549	PT1509228	15	9	22	8
11020550	PT0711228	7	11	22	8
11020551	PT0811228	8	11	22	8
11020552	PT0911228	9	11	22	8
11020553	PT1011228	10	11	22	8
11020554	PT1111228	11	11	22	8
11020555	PT1211228	12	11	22	8
11020556	PT1311228	13	11	22	8
11020557	PT1411228	14	11	22	8
11020558	PT1511228	15	11	22	8




Reference No.	Model No.	H(mm)	W(mm)	L(mm)	Θ°
11020559	PT0709240	7	9	24	0
11020560	PT0809240	8	9	24	0
11020561	PT0909240	9	9	24	0
11020562	PT1009240	10	9	24	0
11020563	PT1109240	11	9	24	0
11020564	PT1209240	12	9	24	0
11020565	PT1309240	13	9	24	0
11020566	PT1409240	14	9	24	0
11020567	PT1509240	15	9	24	0
11020568	PT0711240	7	11	24	0
11020569	PT0811240	8	11	24	0
11020570	PT0911240	9	11	24	0
11020571	PT1011240	10	11	24	0
11020572	PT1111240	11	11	24	0
11020573	PT1211240	12	11	24	0
11020574	PT1311240	13	11	24	0
11020575	PT1411240	14	11	24	0
11020576	PT1511240	15	11	24	0
11020577	PT0709244	7	9	24	4
11020578	PT0809244	8	9	24	4
11020579	PT0909244	9	9	24	4
11020580	PT1009244	10	9	24	4
11020581	PT1109244	11	9	24	4
11020582	PT1209244	12	9	24	4
11020583	PT1309244	13	9	24	4
11020584	PT1409244	14	9	24	4
11020585	PT1509244	15	9	24	4
11020586	PT0711244	7	11	24	4
11020587	PT0811244	8	11	24	4
11020588	PT0911244	9	11	24	4
11020589	PT1011244	10	11	24	4
11020590	PT1111244	11	11	24	4
11020591	PT1211244	12	11	24	4
11020592	PT1311244	13	11	24	4
11020593	PT1411244	14	11	24	4
11020594	PT1511244	15	11	24	4
11020595	PT0709248	7	9	24	8
11020596	PT0809248	8	9	24	8
11020597	PT0909248	9	9	24	8
11020598	PT1009248	10	9	24	8
11020599	PT1109248	11	9	24	8
11020601	PT1209248	12	9	24	8
11020602	PT1309248	13	9	24	8
11020603	PT1409248	14	9	24	8
11020604	PT1509248	15	9	24	8
11020605	PT0711248	7	11	24	8
11020606	PT0811248	8	11	24	8
11020607	PT0911248	9	11	24	8
11020608	PT1011248	10	11	24	8
11020609	PT1111248	11	11	24	8
11020610	PT1211248	12	11	24	8
11020611	PT1311248	13	11	24	8
11020612	PT1411248	14	11	24	8
11020613	PT1511248	15	11	24	8




Reference No.	Model No.	H(mm)	W(mm)	L(mm)	Θ°
11020614	PT0709260	7	9	26	0
11020615	PT0809260	8	9	26	0
11020615	PT0909260	9	9	26	0
11020616	PT1009260	10	9	26	0
11020617	PT1109260	11	9	26	0
11020618	PT1209260	12	9	26	0
11020619	PT1309260	13	9	26	0
11020620	PT1409260	14	9	26	0
11020621	PT1509260	15	9	26	0
11020622	PT0711260	7	11	26	0
11020623	PT0811260	8	11	26	0
11020624	PT0911260	9	11	26	0
11020625	PT1011260	10	11	26	0
11020626	PT1111260	11	11	26	0
11020627	PT1211260	12	11	26	0
11020628	PT1311260	13	11	26	0
11020629	PT1411260	14	11	26	0
11020630	PT1511260	15	11	26	0
11020631	PT0709264	7	9	26	4
11020632	PT0809264	8	9	26	4
11020633	PT0909264	9	9	26	4
11020634	PT1009264	10	9	26	4
11020635	PT1109264	11	9	26	4
11020636	PT1209264	12	9	26	4
11020637	PT1309264	13	9	26	4
11020638	PT1409264	14	9	26	4
11020639	PT1509264	15	9	26	4
11020640	PT0711264	7	11	26	4
11020641	PT0811264	8	11	26	4
11020642	PT0911264	9	11	26	4
11020643	PT1011264	10	11	26	4
11020644	PT1111264	11	11	26	4
11020645	PT1211264	12	11	26	4
11020646	PT1311264	13	11	26	4
11020647	PT1411264	14	11	26	4
11020648	PT1511264	15	11	26	4
11020649	PT0709268	7	9	26	8
11020650	PT0809268	8	9	26	8
11020651	PT0909268	9	9	26	8
11020652	PT1009268	10	9	26	8
11020653	PT1109268	11	9	26	8
11020654	PT1209268	12	9	26	8
11020655	PT1309268	13	9	26	8
11020656	PT1409268	14	9	26	8
11020657	PT1509268	15	9	26	8
11020658	PT0711268	7	11	26	8
11020659	PT0811268	8	11	26	8
11020660	PT0911268	9	11	26	8
11020661	PT1011268	10	11	26	8
11020662	PT1111268	11	11	26	8
11020663	PT1211268	12	11	26	8
11020664	PT1311268	13	11	26	8
11020665	PT1411268	14	11	26	8
11020666	PT1511268	15	11	26	8



Reference No.	Model No.	H(mm)	W(mm)	L(mm)	Θ°
11020667	PT0709280	7	9	28	0
11020668	PT0809280	8	9	28	0
11020669	PT0909280	9	9	28	0
11020670	PT1009280	10	9	28	0
11020671	PT1109280	11	9	28	0
11020672	PT1209280	12	9	28	0
11020673	PT1309280	13	9	28	0
11020674	PT1409280	14	9	28	0
11020675	PT1509280	15	9	28	0
11020676	PT0711280	7	11	28	0
11020677	PT0811280	8	11	28	0
11020678	PT0911280	9	11	28	0
11020679	PT1011280	10	11	28	0
11020680	PT1111280	11	11	28	0
11020681	PT1211280	12	11	28	0
11020682	PT1311280	13	11	28	0
11020683	PT1411280	14	11	28	0
11020684	PT1511280	15	11	28	0
11020685	PT0709284	7	9	28	4
11020686	PT0809284	8	9	28	4
11020687	PT0909284	9	9	28	4
11020688	PT1009284	10	9	28	4
11020689	PT1109284	11	9	28	4
11020690	PT1209284	12	9	28	4
11020691	PT1309284	13	9	28	4
11020692	PT1409284	14	9	28	4
11020693	PT1509284	15	9	28	4
11020694	PT0711284	7	11	28	4
11020695	PT0811284	8	11	28	4
11020696	PT0911284	9	11	28	4
11020697	PT1011284	10	11	28	4
11020698	PT1111284	11	11	28	4
11020699	PT1211284	12	11	28	4
11020700	PT1311284	13	11	28	4
11020701	PT1411284	14	11	28	4
11020702	PT1511284	15	11	28	4
11020703	PT0709288	7	9	28	8
11020704	PT0809288	8	9	28	8
11020705	PT0909288	9	9	28	8
11020706	PT1009288	10	9	28	8
11020707	PT1109288	11	9	28	8
11020708	PT1209288	12	9	28	8
11020709	PT1309288	13	9	28	8
11020710	PT1409288	14	9	28	8
11020711	PT1509288	15	9	28	8
11020712	PT0711288	7	11	28	8
11020713	PT0811288	8	11	28	8
11020714	PT0911288	9	11	28	8
11020715	PT1011288	10	11	28	8
11020716	PT1111288	11	11	28	8
11020717	PT1211288	12	11	28	8
11020718	PT1311288	13	11	28	8
11020719	PT1411288	14	11	28	8
11020720	PT1511288	15	11	28	8




Reference No.	Model No.	H(mm)	W(mm)	L(mm)	Θ°
11020721	PT0709300	7	9	30	0
11020722	PT0809300	8	9	30	0
11020723	PT0909300	9	9	30	0
11020724	PT1009300	10	9	30	0
11020725	PT1109300	11	9	30	0
11020726	PT1209300	12	9	30	0
11020727	PT1309300	13	9	30	0
11020728	PT1409300	14	9	30	0
11020729	PT1509300	15	9	30	0
11020730	PT0711300	7	11	30	0
11020731	PT0811300	8	11	30	0
11020732	PT0911300	9	11	30	0
11020733	PT1011300	10	11	30	0
11020734	PT1111300	11	11	30	0
11020735	PT1211300	12	11	30	0
11020736	PT1311300	13	11	30	0
11020737	PT1411300	14	11	30	0
11020738	PT1511300	15	11	30	0
11020739	PT0709304	7	9	30	4
11020740	PT0809304	8	9	30	4
11020741	PT0909304	9	9	30	4
11020742	PT1009304	10	9	30	4
11020743	PT1109304	11	9	30	4
11020744	PT1209304	12	9	30	4
11020745	PT1309304	13	9	30	4
11020746	PT1409304	14	9	30	4
11020747	PT1509304	15	9	30	4
11020748	PT0711304	7	11	30	4
11020749	PT0811304	8	11	30	4
11020750	PT0911304	9	11	30	4
11020751	PT1011304	10	11	30	4
11020752	PT1111304	11	11	30	4
11020753	PT1211304	12	11	30	4
11020754	PT1311304	13	11	30	4
11020755	PT1411304	14	11	30	4
11020756	PT1511304	15	11	30	4
11020757	PT0709308	7	9	30	8
11020758	PT0809308	8	9	30	8
11020759	PT0909308	9	9	30	8
11020760	PT1009308	10	9	30	8
11020761	PT1109308	11	9	30	8
11020762	PT1209308	12	9	30	8
11020763	PT1309308	13	9	30	8
11020764	PT1409308	14	9	30	8
11020765	PT1509308	15	9	30	8
11020766	PT0711308	7	11	30	8
11020767	PT0811308	8	11	30	8
11020768	PT0911308	9	11	30	8
11020769	PT1011308	10	11	30	8
11020770	PT1111308	11	11	30	8
11020771	PT1211308	12	11	30	8
11020772	PT1311308	13	11	30	8
11020773	PT1411308	14	11	30	8
11020774	PT1511308	15	11	30	8




Reference No.	Model No.	H(mm)	W(mm)	L(mm)	Θ°
11020775	PT0709320	7	9	32	0
11020776	PT0809320	8	9	32	0
11020777	PT0909320	9	9	32	0
11020778	PT1009320	10	9	32	0
11020779	PT1109320	11	9	32	0
11020780	PT1209320	12	9	32	0
11020781	PT1309320	13	9	32	0
11020782	PT1409320	14	9	32	0
11020783	PT1509320	15	9	32	0
11020784	PT0711320	7	11	32	0
11020785	PT0811320	8	11	32	0
11020786	PT0911320	9	11	32	0
11020787	PT1011320	10	11	32	0
11020788	PT1111320	11	11	32	0
11020789	PT1211320	12	11	32	0
11020790	PT1311320	13	11	32	0
11020791	PT1411320	14	11	32	0
11020792	PT1511320	15	11	32	0
11020793	PT0709324	7	9	32	4
11020794	PT0809324	8	9	32	4
11020795	PT0909324	9	9	32	4
11020796	PT1009324	10	9	32	4
11020797	PT1109324	11	9	32	4
11020798	PT1209324	12	9	32	4
11020799	PT1309324	13	9	32	4
11020800	PT1409324	14	9	32	4
11020801	PT1509324	15	9	32	4
11020802	PT0711324	7	11	32	4
11020803	PT0811324	8	11	32	4
11020804	PT0911324	9	11	32	4
11020805	PT1011324	10	11	32	4
11020806	PT1111324	11	11	32	4
11020807	PT1211324	12	11	32	4
11020808	PT1311324	13	11	32	4
11020809	PT1411324	14	11	32	4
11020810	PT1511324	15	11	32	4
11020811	PT0709328	7	9	32	8
11020812	PT0809328	8	9	32	8
11020813	PT0909328	9	9	32	8
11020814	PT1009328	10	9	32	8
11020815	PT1109328	11	9	32	8
11020816	PT1209328	12	9	32	8
11020817	PT1309328	13	9	32	8
11020818	PT1409328	14	9	32	8
11020819	PT1509328	15	9	32	8
11020820	PT0711328	7	11	32	8
11020821	PT0811328	8	11	32	8
11020822	PT0911328	9	11	32	8
11020823	PT1011328	10	11	32	8
11020824	PT1111328	11	11	32	8
11020825	PT1211328	12	11	32	8
11020826	PT1311328	13	11	32	8
11020827	PT1411328	14	11	32	8
11020828	PT1511328	15	11	32	8


3d Cage-TLIF sterile packaged cages




Reference No.	Model No.	H(mm)	W(mm)	L(mm)	Θ°
11020829	TT0710280	7	10	28	0
11020830	TT0810280	8	10	28	0
11020831	TT0910280	9	10	28	0
11020832	TT1010280	10	10	28	0
11020833	TT1110280	11	10	28	0
11020834	TT1210280	12	10	28	0
11020835	TT1310280	13	10	28	0
11020836	TT1410280	14	10	28	0
11020837	TT1510280	15	10	28	0
11020838	TT0712280	7	12	28	0
11020839	TT0812280	8	12	28	0
11020840	TT0912280	9	12	28	0
11020841	TT1012280	10	12	28	0
11020842	TT1112280	11	12	28	0
11020843	TT1212280	12	12	28	0
11020844	TT1312280	13	12	28	0
11020845	TT1412280	14	12	28	0
11020846	TT1512280	15	12	28	0
11020847	TT0710284	7	10	28	4
11020848	TT0810284	8	10	28	4
11020849	TT0910284	9	10	28	4
11020850	TT1010284	10	10	28	4
11020851	TT1110284	11	10	28	4
11020852	TT1210284	12	10	28	4
11020853	TT1310284	13	10	28	4
11020854	TT1410284	14	10	28	4
11020855	TT1510284	15	10	28	4
11020856	TT0712284	7	12	28	4
11020857	TT0812284	8	12	28	4
11020858	TT0912284	9	12	28	4
11020859	TT1012284	10	12	28	4
11020860	TT1112284	11	12	28	4
11020861	TT1212284	12	12	28	4
11020862	TT1312284	13	12	28	4
11020863	TT1412284	14	12	28	4
11020864	TT1512284	15	12	28	4
11020865	TT0710288	7	10	28	8
11020866	TT0810288	8	10	28	8
11020867	TT0910288	9	10	28	8
11020868	TT1010288	10	10	28	8
11020869	TT1110288	11	10	28	8
11020870	TT1210288	12	10	28	8
11020871	TT1310288	13	10	28	8
11020872	TT1410288	14	10	28	8
11020873	TT1510288	15	10	28	8
11020874	TT0712288	7	12	28	8
11020875	TT0812288	8	12	28	8
11020876	TT0912288	9	12	28	8
11020877	TT1012288	10	12	28	8
11020878	TT1112288	11	12	28	8
11020879	TT1212288	12	12	28	8
11020880	TT1312288	13	12	28	8
11020881	TT1412288	14	12	28	8
11020882	TT1512288	15	12	28	8
11020883	TT0710282	7	10	28	12
11020884	TT0810282	8	10	28	12
11020885	TT0910282	9	10	28	12
11020886	TT1010282	10	10	28	12
11020887	TT1110282	11	10	28	12
11020888	TT1210282	12	10	28	12
11020889	TT1310282	13	10	28	12
11020890	TT1410282	14	10	28	12
11020891	TT1510282	15	10	28	12
11020892	TT0712282	7	12	28	12
11020893	TT0812282	8	12	28	12
11020894	TT0912282	9	12	28	12
11020895	TT1012282	10	12	28	12
11020896	TT1112282	11	12	28	12
11020897	TT1212282	12	12	28	12
11020898	TT1312282	13	12	28	12
11020899	TT1412282	14	12	28	12
11020900	TT1512282	15	12	28	12




Reference No.	Model No.	H(mm)	W(mm)	L(mm)	Θ°
11020901	TT0710300	7	10	30	0
11020902	TT0810300	8	10	30	0
11020903	TT0910300	9	10	30	0
11020904	TT1010300	10	10	30	0
11020905	TT1110300	11	10	30	0
11020906	TT1210300	12	10	30	0
11020907	TT1310300	13	10	30	0
11020908	TT1410300	14	10	30	0
11020909	TT1510300	15	10	30	0
11020910	TT0712300	7	12	30	0
11020911	TT0812300	8	12	30	0
11020912	TT0912300	9	12	30	0
11020913	TT1012300	10	12	30	0
11020914	TT1112300	11	12	30	0
11020915	TT1212300	12	12	30	0
11020916	TT1312300	13	12	30	0
11020917	TT1412300	14	12	30	0
11020918	TT1512300	15	12	30	0
11020919	TT0710304	7	10	30	4
11020920	TT0810304	8	10	30	4
11020921	TT0910304	9	10	30	4
11020922	TT1010304	10	10	30	4
11020923	TT1110304	11	10	30	4
11020924	TT1210304	12	10	30	4
11020925	TT1310304	13	10	30	4
11020926	TT1410304	14	10	30	4
11020927	TT1510304	15	10	30	4
11020928	TT0712304	7	12	30	4
11020929	TT0812304	8	12	30	4
11020930	TT0912304	9	12	30	4
11020931	TT1012304	10	12	30	4
11020932	TT1112304	11	12	30	4
11020933	TT1212304	12	12	30	4
11020934	TT1312304	13	12	30	4
11020935	TT1412304	14	12	30	4
11020936	TT1512304	15	12	30	4
11020937	TT0710308	7	10	30	8
11020938	TT0810308	8	10	30	8
11020939	TT0910308	9	10	30	8
11020940	TT1010308	10	10	30	8
11020941	TT1110308	11	10	30	8
11020942	TT1210308	12	10	30	8
11020943	TT1310308	13	10	30	8
11020944	TT1410308	14	10	30	8
11020945	TT1510308	15	10	30	8
11020946	TT0712308	7	12	30	8
11020947	TT0812308	8	12	30	8
11020948	TT0912308	9	12	30	8
11020949	TT1012308	10	12	30	8
11020950	TT1112308	11	12	30	8
11020951	TT1212308	12	12	30	8
11020952	TT1312308	13	12	30	8
11020953	TT1412308	14	12	30	8
11020954	TT1512308	15	12	30	8
11020955	TT0710302	7	10	30	12
11020956	TT0810302	8	10	30	12
11020957	TT0910302	9	10	30	12
11020958	TT1010302	10	10	30	12
11020959	TT1110302	11	10	30	12
11020960	TT1210302	12	10	30	12
11020961	TT1310302	13	10	30	12
11020962	TT1410302	14	10	30	12
11020963	TT1510302	15	10	30	12
11020964	TT0712302	7	12	30	12
11020965	TT0812302	8	12	30	12
11020966	TT0912302	9	12	30	12
11020967	TT1012302	10	12	30	12
11020968	TT1112302	11	12	30	12
11020969	TT1212302	12	12	30	12
11020970	TT1312302	13	12	30	12
11020971	TT1412302	14	12	30	12
11020972	TT1512302	15	12	30	12




Reference No.	Model No.	H(mm)	W(mm)	L(mm)	Θ°
11020971	TT0710320	7	10	32	0
11020972	TT0810320	8	10	32	0
11020973	TT0910320	9	10	32	0
11020974	TT1010320	10	10	32	0
11020975	TT1110320	11	10	32	0
11020976	TT1210320	12	10	32	0
11020977	TT1310320	13	10	32	0
11020978	TT1410320	14	10	32	0
11020979	TT1510320	15	10	32	0
11020980	TT0712320	7	12	32	0
11020981	TT0812320	8	12	32	0
11020982	TT0912320	9	12	32	0
11020983	TT1012320	10	12	32	0
11020984	TT1112320	11	12	32	0
11020985	TT1212320	12	12	32	0
11020986	TT1312320	13	12	32	0
11020987	TT1412320	14	12	32	0
11020988	TT1512320	15	12	32	0
11020989	TT0710324	7	10	32	4
11020990	TT0810324	8	10	32	4
11020991	TT0910324	9	10	32	4
11020992	TT1010324	10	10	32	4
11020993	TT1110324	11	10	32	4
11020994	TT1210324	12	10	32	4
11020995	TT1310324	13	10	32	4
11020996	TT1410324	14	10	32	4
11020997	TT1510324	15	10	32	4
11020998	TT0712324	7	12	32	4
11020999	TT0812324	8	12	32	4
11021001	TT0912324	9	12	32	4
11021002	TT1012324	10	12	32	4
11021003	TT1112324	11	12	32	4
11021004	TT1212324	12	12	32	4
11021005	TT1312324	13	12	32	4
11021006	TT1412324	14	12	32	4
11021007	TT1512324	15	12	32	4
11021008	TT0710328	7	10	32	8
11021009	TT0810328	8	10	32	8
11021010	TT0910328	9	10	32	8
11021011	TT1010328	10	10	32	8
11021012	TT1110328	11	10	32	8
11021013	TT1210328	12	10	32	8
11021014	TT1310328	13	10	32	8
11021015	TT1410328	14	10	32	8
11021016	TT1510328	15	10	32	8
11021017	TT0712328	7	12	32	8
11021018	TT0812328	8	12	32	8
11021019	TT0912328	9	12	32	8
11021020	TT1012328	10	12	32	8
11021021	TT1112328	11	12	32	8
11021022	TT1212328	12	12	32	8
11021023	TT1312328	13	12	32	8
11021024	TT1412328	14	12	32	8
11021025	TT1512328	15	12	32	8
11021026	TT0710322	7	10	32	12
11021027	TT0810322	8	10	32	12
11021028	TT0910322	9	10	32	12
11021029	TT1010322	10	10	32	12
11021030	TT1110322	11	10	32	12
11021031	TT1210322	12	10	32	12
11021032	TT1310322	13	10	32	12
11021033	TT1410322	14	10	32	12
11021034	TT1510322	15	10	32	12
11021035	TT0712322	7	12	32	12
11021036	TT0812322	8	12	32	12
11021037	TT0912322	9	12	32	12
11021038	TT1012322	10	12	32	12
11021039	TT1112322	11	12	32	12
11021040	TT1212322	12	12	32	12
11021041	TT1312322	13	12	32	12
11021042	TT1412322	14	12	32	12
11021043	TT1512322	15	12	32	12



Reference No.	Model No.	H(mm)	W(mm)	L(mm)	Θ°
11021044	TT0710340	7	10	34	0
11021045	TT0810340	8	10	34	0
11021046	TT0910340	9	10	34	0
11021047	TT1010340	10	10	34	0
11021048	TT1110340	11	10	34	0
11021049	TT1210340	12	10	34	0
11021050	TT1310340	13	10	34	0
11021051	TT1410340	14	10	34	0
11021052	TT1510340	15	10	34	0
11021053	TT0712340	7	12	34	0
11021054	TT0812340	8	12	34	0
11021055	TT0912340	9	12	34	0
11021056	TT1012340	10	12	34	0
11021057	TT1112340	11	12	34	0
11021058	TT1212340	12	12	34	0
11021059	TT1312340	13	12	34	0
11021060	TT1412340	14	12	34	0
11021061	TT1512340	15	12	34	0
11021062	TT0710344	7	10	34	4
11021063	TT0810344	8	10	34	4
11021064	TT0910344	9	10	34	4
11021065	TT1010344	10	10	34	4
11021066	TT1110344	11	10	34	4
11021067	TT1210344	12	10	34	4
11021068	TT1310344	13	10	34	4
11021069	TT1410344	14	10	34	4
11021070	TT1510344	15	10	34	4
11021071	TT0712344	7	12	34	4
11021072	TT0812344	8	12	34	4
11021073	TT0912344	9	12	34	4
11021074	TT1012344	10	12	34	4
11021075	TT1112344	11	12	34	4
11021076	TT1212344	12	12	34	4
11021077	TT1312344	13	12	34	4
11021078	TT1412344	14	12	34	4
11021079	TT1512344	15	12	34	4
11021080	TT0710348	7	10	34	8
11021081	TT0810348	8	10	34	8
11021082	TT0910348	9	10	34	8
11021083	TT1010348	10	10	34	8
11021084	TT1110348	11	10	34	8
11021085	TT1210348	12	10	34	8
11021086	TT1310348	13	10	34	8
11021087	TT1410348	14	10	34	8
11021088	TT1510348	15	10	34	8
11021089	TT0712348	7	12	34	8
11021090	TT0812348	8	12	34	8
11021091	TT0912348	9	12	34	8
11021092	TT1012348	10	12	34	8
11021093	TT1112348	11	12	34	8
11021094	TT1212348	12	12	34	8
11021095	TT1312348	13	12	34	8
11021096	TT1412348	14	12	34	8
11021097	TT1512348	15	12	34	8
11021098	TT0710342	7	10	34	12
11021099	TT0810342	8	10	34	12
11021101	TT0910342	9	10	34	12
11021102	TT1010342	10	10	34	12
11021103	TT1110342	11	10	34	12
11021104	TT1210342	12	10	34	12
11021105	TT1310342	13	10	34	12
11021106	TT1410342	14	10	34	12
11021107	TT1510342	15	10	34	12
11021108	TT0712342	7	12	34	12
11021109	TT0812342	8	12	34	12
11021110	TT0912342	9	12	34	12
11021111	TT1012342	10	12	34	12
11021112	TT1112342	11	12	34	12
11021113	TT1212342	12	12	34	12
11021114	TT1312342	13	12	34	12
11021115	TT1412342	14	12	34	12
11021116	TT1512342	15	12	34	12




Reference No.	Model No.	H(mm)	W(mm)	L(mm)	Θ°
11021117	TT0710360	7	10	36	0
11021118	TT0810360	8	10	36	0
11021119	TT0910360	9	10	36	0
11021120	TT1010360	10	10	36	0
11021121	TT1110360	11	10	36	0
11021122	TT1210360	12	10	36	0
11021123	TT1310360	13	10	36	0
11021124	TT1410360	14	10	36	0
11021125	TT1510360	15	10	36	0
11021126	TT0712360	7	12	36	0
11021127	TT0812360	8	12	36	0
11021128	TT0912360	9	12	36	0
11021129	TT1012360	10	12	36	0
11021130	TT1112360	11	12	36	0
11021131	TT1212360	12	12	36	0
11021132	TT1312360	13	12	36	0
11021133	TT1412360	14	12	36	0
11021134	TT1512360	15	12	36	0
11021135	TT0710364	7	10	36	4
11021136	TT0810364	8	10	36	4
11021137	TT0910364	9	10	36	4
11021138	TT1010364	10	10	36	4
11021139	TT1110364	11	10	36	4
11021140	TT1210364	12	10	36	4
11021141	TT1310364	13	10	36	4
11021142	TT1410364	14	10	36	4
11021143	TT1510364	15	10	36	4
11021144	TT0712364	7	12	36	4
11021145	TT0812364	8	12	36	4
11021146	TT0912364	9	12	36	4
11021147	TT1012364	10	12	36	4
11021148	TT1112364	11	12	36	4
11021149	TT1212364	12	12	36	4
11021150	TT1312364	13	12	36	4
11021151	TT1412364	14	12	36	4
11021152	TT1512364	15	12	36	4
11021153	TT0710368	7	10	36	8
11021154	TT0810368	8	10	36	8
11021155	TT0910368	9	10	36	8
11021156	TT1010368	10	10	36	8
11021157	TT1110368	11	10	36	8
11021158	TT1210368	12	10	36	8
11021159	TT1310368	13	10	36	8
11021160	TT1410368	14	10	36	8
11021161	TT1510368	15	10	36	8
11021162	TT0712368	7	12	36	8
11021163	TT0812368	8	12	36	8
11021164	TT0912368	9	12	36	8
11021165	TT1012368	10	12	36	8
11021166	TT1112368	11	12	36	8
11021167	TT1212368	12	12	36	8
11021168	TT1312368	13	12	36	8
11021169	TT1412368	14	12	36	8
11021170	TT1512368	15	12	36	8
11021171	TT0710362	7	10	36	12
11021172	TT0810362	8	10	36	12
11021173	TT0910362	9	10	36	12
11021174	TT1010362	10	10	36	12
11021175	TT1110362	11	10	36	12
11021176	TT1210362	12	10	36	12
11021177	TT1310362	13	10	36	12
11021178	TT1410362	14	10	36	12
11021179	TT1510362	15	10	36	12
11021180	TT0712362	7	12	36	12
11021181	TT0812362	8	12	36	12
11021182	TT0912362	9	12	36	12
11021183	TT1012362	10	12	36	12
11021184	TT1112362	11	12	36	12
11021185	TT1212362	12	12	36	12
11021186	TT1312362	13	12	36	12
11021187	TT1412362	14	12	36	12
11021188	TT1512362	15	12	36	12



Reference No.

11021189
11021190
11021191
11021192
11021193
11021194
11021195
11021196
11021197
11021198
11021199
11021201
11021202
11021203
11021204
11021205
11021206
11021207
11021208
11021209
11021210
11021211
11021212
11021213
11021214
11021215
11021216
11021217
11021218
11021219
11021220
11021221
11021222
11021223
11021224
11021225
11021226
11021227
11021228
11021229
11021230
11021231
11021232
11021233
11021234
11021235
11021236
11021237
11021238
11021239
11021240
11021241
11021242
11021243
11021244
11021245
11021246
11021247
11021248
11021249
11021250
11021251
11021252
11021253
11021254
11021255
11021256
11021257
11021258
11021259
11021260
11021261

Model No.	H(mm)	W(mm)	L(mm)	Θ°
TT0710380	7	10	38	0
TT0810380	8	10	38	0
TT0910380	9	10	38	0
TT1010380	10	10	38	0
TT1110380	11	10	38	0
TT1210380	12	10	38	0
TT1310380	13	10	38	0
TT1410380	14	10	38	0
TT1510380	15	10	38	0
TT0712380	7	12	38	0
TT0812380	8	12	38	0
TT0912380	9	12	38	0
TT1012380	10	12	38	0
TT1112380	11	12	38	0
TT1212380	12	12	38	0
TT1312380	13	12	38	0
TT1412380	14	12	38	0
TT1512380	15	12	38	0
TT0710384	7	10	38	4
TT0810384	8	10	38	4
TT0910384	9	10	38	4
TT1010384	10	10	38	4
TT1110384	11	10	38	4
TT1210384	12	10	38	4
TT1310384	13	10	38	4
TT1410384	14	10	38	4
TT1510384	15	10	38	4
TT0712384	7	12	38	4
TT0812384	8	12	38	4
TT0912384	9	12	38	4
TT1012384	10	12	38	4
TT1112384	11	12	38	4
TT1212384	12	12	38	4
TT1312384	13	12	38	4
TT1412384	14	12	38	4
TT1512384	15	12	38	4
TT0710388	7	10	38	8
TT0810388	8	10	38	8
TT0910388	9	10	38	8
TT1010388	10	10	38	8
TT1110388	11	10	38	8
TT1210388	12	10	38	8
TT1310388	13	10	38	8
TT1410388	14	10	38	8
TT1510388	15	10	38	8
TT0712388	7	12	38	8
TT0812388	8	12	38	8
TT0912388	9	12	38	8
TT1012388	10	12	38	8
TT1112388	11	12	38	8
TT1212388	12	12	38	8
TT1312388	13	12	38	8
TT1412388	14	12	38	8
TT1512388	15	12	38	8
TT0710382	7	10	38	12
TT0810382	8	10	38	12
TT0910382	9	10	38	12
TT1010382	10	10	38	12
TT1110382	11	10	38	12
TT1210382	12	10	38	12
TT1310382	13	10	38	12
TT1410382	14	10	38	12
TT1510382	15	10	38	12
TT0712382	7	12	38	12
TT0812382	8	12	38	12
TT0912382	9	12	38	12
TT1012382	10	12	38	12
TT1112382	11	12	38	12
TT1212382	12	12	38	12
TT1312382	13	12	38	12
TT1412382	14	12	38	12
TT1512382	15	12	38	12




Reference No.	Model No.	H(mm)	W(mm)	L(mm)	Θ°
11021262	TT0710400	7	10	40	0
11021263	TT0810400	8	10	40	0
11021264	TT0910400	9	10	40	0
11021265	TT1010400	10	10	40	0
11021266	TT1110400	11	10	40	0
11021267	TT1210400	12	10	40	0
11021268	TT1310400	13	10	40	0
11021269	TT1410400	14	10	40	0
11021270	TT1510400	15	10	40	0
11021271	TT0712400	7	12	40	0
11021272	TT0812400	8	12	40	0
11021273	TT0912400	9	12	40	0
11021274	TT1012400	10	12	40	0
11021275	TT1112400	11	12	40	0
11021276	TT1212400	12	12	40	0
11021277	TT1312400	13	12	40	0
11021278	TT1412400	14	12	40	0
11021279	TT1512400	15	12	40	0
11021280	TT0710404	7	10	40	4
11021281	TT0810404	8	10	40	4
11021282	TT0910404	9	10	40	4
11021283	TT1010404	10	10	40	4
11021284	TT1110404	11	10	40	4
11021285	TT1210404	12	10	40	4
11021286	TT1310404	13	10	40	4
11021287	TT1410404	14	10	40	4
11021288	TT1510404	15	10	40	4
11021289	TT0712404	7	12	40	4
11021290	TT0812404	8	12	40	4
11021291	TT0912404	9	12	40	4
11021292	TT1012404	10	12	40	4
11021293	TT1112404	11	12	40	4
11021294	TT1212404	12	12	40	4
11021295	TT1312404	13	12	40	4
11021296	TT1412404	14	12	40	4
11021297	TT1512404	15	12	40	4
11021298	TT0710408	7	10	40	8
11021299	TT0810408	8	10	40	8
11021301	TT0910408	9	10	40	8
11021302	TT1010408	10	10	40	8
11021303	TT1110408	11	10	40	8
11021304	TT1210408	12	10	40	8
11021305	TT1310408	13	10	40	8
11021306	TT1410408	14	10	40	8
11021307	TT1510408	15	10	40	8
11021308	TT0712408	7	12	40	8
11021309	TT0812408	8	12	40	8
11021310	TT0912408	9	12	40	8
11021311	TT1012408	10	12	40	8
11021312	TT1112408	11	12	40	8
11021313	TT1212408	12	12	40	8
11021314	TT1312408	13	12	40	8
11021315	TT1412408	14	12	40	8
11021316	TT1512408	15	12	40	8
11021317	TT0710402	7	10	40	12
11021318	TT0810402	8	10	40	12
11021319	TT0910402	9	10	40	12
11021320	TT1010402	10	10	40	12
11021311	TT1110402	11	10	40	12
11021312	TT1210402	12	10	40	12
11021313	TT1310402	13	10	40	12
11021314	TT1410402	14	10	40	12
11021315	TT1510402	15	10	40	12
11021316	TT0712402	7	12	40	12
11021317	TT0812402	8	12	40	12
11021318	TT0912402	9	12	40	12
11021319	TT1012402	10	12	40	12
11021320	TT1112402	11	12	40	12
11021321	TT1212402	12	12	40	12
11021322	TT1312402	13	12	40	12
11021323	TT1412402	14	12	40	12
11021324	TT1512402	15	12	40	12

GENOSS

Lumbar Pedicle Screw

Order Information

Lumbar Pedicle Screw sterile & non-sterile packaged screws



W : S.L.A. Surface treatment (Sterilized)
N : Non S.L.A. Surface treatment (Non-sterile)

No.	REF	D(Ø)	L(mm)
1	LSP4025SW, N	4.0	25
2	LSP4030SW, N	4.0	30
3	LSP4035SW, N	4.0	35
4	LSP4040SW, N	4.0	40
5	LSP4045SW, N	4.0	45
6	LSP4050SW, N	4.0	50
7	LSP4530SW, N	4.5	30
8	LSP4535SW, N	4.5	35
9	LSP4540SW, N	4.5	40
10	LSP4545SW, N	4.5	45
11	LSP4550SW, N	4.5	50
12	LSP4555SW, N	4.5	55
13	LSP4560SW, N	4.5	60
14	LSP4565SW, N	4.5	65
15	LSP5030SW, N	5.0	30
16	LSP5035SW, N	5.0	35
17	LSP5040SW, N	5.0	40
18	LSP5045SW, N	5.0	45
19	LSP5050SW, N	5.0	50
20	LSP5055SW, N	5.0	55
21	LSP5060SW, N	5.0	60
22	LSP5065SW, N	5.0	65
23	LSP5530SW, N	5.5	30
24	LSP5535SW, N	5.5	35

No.	REF	D(Ø)	L(mm)
1	LSP5540SW, N	5.5	40
2	LSP5545SW, N	5.5	45
3	LSP5550SW, N	5.5	50
4	LSP5555SW, N	5.5	55
5	LSP5560SW, N	5.5	60
6	LSP5565SW, N	5.5	65
7	LSP6030SW, N	6.0	30
8	LSP6035SW, N	6.0	35
9	LSP6040SW, N	6.0	40
10	LSP6045SW, N	6.0	45
11	LSP6050SW, N	6.0	50
12	LSP6055SW, N	6.0	55
13	LSP6060SW, N	6.0	60
14	LSP6065SW, N	6.0	65
15	LSP6530SW, N	6.5	30
16	LSP6535SW, N	6.5	35
17	LSP6540SW, N	6.5	40
18	LSP6545SW, N	6.5	45
19	LSP6550SW, N	6.5	50
20	LSP6555SW, N	6.5	55
21	LSP6560SW, N	6.5	60
22	LSP6565SW, N	6.5	65
23	LSP7030SW, N	7.0	30
24	LSP7035SW, N	7.0	35

No.	REF	D(Ø)	L(mm)
1	LSP7040SW, N	7.0	40
2	LSP7045SW, N	7.0	45
3	LSP7050SW, N	7.0	50
4	LSP7055SW, N	7.0	55
5	LSP7060SW, N	7.0	60
6	LSP7065SW, N	7.0	65
7	LSP7530SW, N	7.5	30
8	LSP7535SW, N	7.5	35
9	LSP7540SW, N	7.5	40
10	LSP7545SW, N	7.5	45
11	LSP7550SW, N	7.5	50
12	LSP7555SW, N	7.5	55
13	LSP7560SW, N	7.5	60
14	LSP7565SW, N	7.5	65
15	LSP8030SW, N	8.0	30
16	LSP8035SW, N	8.0	35
17	LSP8040SW, N	8.0	40
18	LSP8045SW, N	8.0	45
19	LSP8050SW, N	8.0	50
20	LSP8055SW, N	8.0	55
21	LSP8060SW, N	8.0	60
22	LSP8065SW, N	8.0	65
23	LSP8070SW, N	8.0	70
24	LSP8075SW, N	8.0	75

GENOSS
Spinal Devices



GENOSS 3d Cage™

Manufacturer FIMS Co., Ltd.

56, Changnyong-daero 256beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do 16229, REPUBLIC OF KOREA

Tel. +82-70-7098-8806 Fax. +82-31-261-8649

Distributor GENOSS Co., Ltd.

1F, Gyeonggi R&DB Center/226, 2F, GSBC, 105 Gwanggyo-ro, Yeongtong-gu, Suwon-si, Gyeonggi-do, Korea
www.genoss.co.kr

Tel. +82-31-888-5100 Fax. +82-31-888-5105

GENOSS Lumbar Pedicle Screw™

Manufacturer GENOSS Co., Ltd.

1F, Gyeonggi R&DB Center/226, 2F, GSBC, 105 Gwanggyo-ro, Yeongtong-gu, Suwon-si, Gyeonggi-do, Korea
www.genoss.co.kr

Tel. +82-31-888-5100 Fax. +82-31-888-5105